


# A first-of-its-kind multi-omic assay using lipids and proteins accurately identifies early-stage ovarian cancer in symptomatic individuals



Check out our poster from AACR 2026 to see how we got here!

Rachel Culp-Hill<sup>1</sup>, Brendan Giles<sup>1</sup>, Shannon Kilkenny<sup>1</sup>, Charles Nichols<sup>1</sup>, Yu Han<sup>1</sup>, Moisés Zapata<sup>1</sup>, Enkhtuya Radnaa<sup>1</sup>, Maria Wong<sup>1</sup>, Mattie Goldberg<sup>1</sup>, Vuna S. Fa<sup>1</sup>, Kim Ekroos<sup>2</sup>, James Robert White<sup>3</sup>, Abigail McElhinny<sup>1</sup>  
<sup>1</sup>AOA Dx, Denver, CO, USA | <sup>2</sup>Lipidomics Consulting, LLC, Espoo, Finland | <sup>3</sup>Resphera Biosciences, Baltimore, MD, USA


Check out our publication in *Cancer Research Communications!*




## Background and Rationale

- Ovarian cancer (OC) remains highly lethal** primarily due to late-stage diagnosis<sup>1</sup>. Existing diagnostic tools lack sensitivity and specificity to detect early-stage disease<sup>2</sup>.
- While stage I disease shows survival rates above 90%, most OC is detected **at late stages (III/IV)**, when **5-year survival is <30%**<sup>1</sup>.
- ~85% of patients present with symptoms<sup>3</sup>**, but they are often vague and misattributed, leading to a 9-month diagnostic delay<sup>4</sup>.
- HGSOC **doubles in volume every 2-4 months**, leaving a narrow detection window before metastasis; earlier diagnosis is associated with significantly improved survival<sup>5,6,7</sup>.
- There is a clear **unmet need for improved diagnostic tools** to accurately triage symptomatic women and enable earlier detection.


## Biomarker Down-selection



## Why Lipids? Biologically Relevant



## Patient Journey to OC Diagnosis



## Study Design + Methods

Group Description	Total
All OC	116
Early-Stage OC (I/II)	50
Late-Stage OC (III/IV)	66
All Controls	230
Benign gyn. conditions + masses	116
Symptomatic	114
<b>Grand Total</b>	<b>346</b>

- Targeted lipid + protein biomarkers evaluated from prior discovery<sup>8</sup>
- Prospective, multi-site U.S. cohort of symptomatic women
- Includes OC + controls (benign + longitudinally confirmed non-malignant)
- Treatment-naïve, blinded samples
- Supplemental samples from commercial source
- Integrated lipid + protein dataset analyzed using multivariable modeling


## Why MS? Proven + Scalable

Mass spectrometry is already widely deployed in clinical labs<sup>10,11</sup>:

Endocrinology	Newborn Screening	Vitamins	Clinical Toxicology
<ul style="list-style-type: none"> <li><b>Comprehensive hormone panels</b> (congenital adrenal hyperplasia)</li> <li><b>Cortisol</b> (Cushing's syndrome)</li> <li><b>Sex steroid panels</b> (hypogonadism, PCOS)</li> </ul>	<ul style="list-style-type: none"> <li><b>Newborn screening panels</b> (&gt;30-60 inborn metabolic errors)</li> <li><b>Acylcarnitine profiling</b> (fatty acid oxidation disorders)</li> <li><b>Amino acid analysis</b> (phenylketonuria)</li> </ul>	<ul style="list-style-type: none"> <li><b>Vitamin D</b> (deficiency diagnosis)</li> <li><b>Vitamin B6</b> (nutritional, metabolic assessment)</li> <li><b>Vitamin A + E</b> (malabsorption, deficiency disorders)</li> </ul>	<ul style="list-style-type: none"> <li><b>Opioid + drug panels</b> (pain management)</li> <li><b>Therapeutic drug monitoring</b> (anti-epileptics)</li> <li><b>Toxicology</b> (overdose evaluation)</li> </ul>

Compatible with standard clinical workflows and automation, providing robust and reproducible measurements.

## Targeted Assay Development



## Clinical Performance in Intended-Use Population

A targeted multi-omic model (**8 lipids + 3 proteins**) demonstrates strong and consistent discrimination of ovarian cancer, including early-stage, in a symptomatic population.

Analyte	Within-Day %CV	Between-Day %CV
Lipid #1	1.84%	0.59%
Lipid #2	0.79%	1.29%
Lipid #3	1.22%	2.43%
Lipid #4	1.59%	3.70%
Lipid #5	1.50%	2.64%
Lipid #6	1.37%	1.10%
Lipid #7	3.05%	3.42%
Lipid #8	4.80%	0.00%

**OC vs. Symptomatic**  
Sensitivity: 96.5%  
Specificity: 80.6%

**Early-stage OC vs. Symptomatic**  
Sensitivity: 92%  
Specificity: 80.6%

Performance is maintained when translated to a targeted, clinical-grade assay, further supporting the performance from previous discovery cohorts (Giles & Culp-Hill et al. 2025)

## Analytical Robustness

The targeted lipid panel demonstrates high analytical precision (CV <5%) across runs, supporting reproducible measurement in a clinical workflow.

Panel features were selected in part for consistent, reliable quantification using targeted LC-MS.

## Conclusions

- High sensitivity detection of early-stage OC** in symptomatic individuals
- Final panel defined: 8 lipids + 3 proteins**, translated into a targeted clinical-grade assay
- Performance holds** in a prospective, real-world intended-use population
- Establishes a clinically aligned, deployable multi-omic diagnostic**

References: 1. Siegel RL, Miller KD, Fuchs HE, et al. Cancer Statistics, 2023. *CA Cancer Clin Oncol*. 2023;73(1):17-48. 2. Henson DE, et al. The Role of the Gb3-CD in Early Detection of Epithelial OC in Women at Average Risk. *Obstet Gynecol*. 2017;129(4):646-53. 3. Cote AG, et al. High-grade serous ovarian cancer: a review. *JAMA Network Open*. 2020;3(10):e201910. 4. Giles B, et al. Early-stage ovarian cancer: a review. *JAMA Network Open*. 2020;3(10):e201910. 5. Giles B, et al. The potential natural history of serous ovarian cancer: defining the target for early detection. *PLoS Med*. 2019;16(10):e1002811. 6. Siegel RL, et al. Cancer Statistics, 2023. *CA Cancer Clin Oncol*. 2023;73(1):17-48. 7. Siegel RL, et al. Cancer Statistics, 2023. *CA Cancer Clin Oncol*. 2023;73(1):17-48. 8. Giles B, et al. Early-stage ovarian cancer: a review. *JAMA Network Open*. 2020;3(10):e201910. 9. Siegel RL, et al. Cancer Statistics, 2023. *CA Cancer Clin Oncol*. 2023;73(1):17-48. 10. Siegel RL, et al. Cancer Statistics, 2023. *CA Cancer Clin Oncol*. 2023;73(1):17-48. 11. Siegel RL, et al. Cancer Statistics, 2023. *CA Cancer Clin Oncol*. 2023;73(1):17-48.